

➤



Stock **#37-330** [CONTACT](#)

1

**€1.412,00**

AJOUTER AU PANIER

**€1.412,00**

## Demande de Devis

! Les prix sont indiqués hors TVA et droits applicables.

Espace téléchargement

**TECHSPEC® 90° Kinematic Adapter (2), TECHSPEC® Coherent® OBIS™ Heat Sink (2), TECHSPEC® 25/25.4mm Diameter Kinematic Mount, 3-Screws (2), TECHSPEC® 30mm ID Standard Cage Plate with 1/4-20 for Post Mounting (6), TECHSPEC® 6mm Diameter x 75mm Length, M3, Cage Support Rod (8), TECHSPEC® 25mm ID Standard Cage Plate with 1/4-20 for Post Mounting (1), TECHSPEC® Cage Sphere with (5) 30mm Ports and (1) 43mm Port (1), TECHSPEC® 25/25.4mm Cube Beamsplitter Mount (1), TECHSPEC® 2" Post Holder and #03-655, Post Holder Assembly (2), TECHSPEC® 3" Length, 1/4-20 Stud, Steel Post (2)**

---

Conformité réglementaire

RoHS:  
[Conforme](#)

Certificate of Conformance:  
[Visionner](#)

---