

[Afficher tous les 17 produits de la même famille.](#)

Olympus SZ51/SZ61 Objectif 0,4X

See More by [Olympus](#)



SZ61 Stereo Microscope (eyepiece, objective, and stand each sold separately)



Stock **#88-128** [CONTACT](#)

-

1

+

€422⁰⁰

AJOUTER AU PANIER

Prix sur Quantité	
Qté 1+	€422,00 prix unitaire
Need More?	Demande de Devis

Les prix sont indiqués hors TVA et droits applicables.

Espace téléchargement

Caractéristiques du produit	
110ALK0.4X/ 1-S818	Numéro de Modèle:
Olympus	Fabricant:
Propriétés optiques	

0.4X

Certificate of Conformance:

Description produit

- Les Microscopes Zoom Stéréo d'Olympus SZ51 et SZ61 sont utilisés dans les sciences de la vie pour offrir une qualité d'image supérieure et proposent une grande profondeur de champ idéale pour de nombreuses applications en microscopie biologique. Ces Microscopes utilisent un système optique Greenough disposant d'un angle de convergence de 10° qui présente une excellente planéité de champ, ainsi qu'une grande clarté, des détails très marqués et des couleurs très précises. Chaque microscope est doté de commandes d'accès frontal pratiques, d'une bague de butée pour le zoom permettant d'effectuer des changements de grossissement rapidement reproductibles et d'un oculaire confortable, disposant d'un contrôle d'aberration de la pupille et d'un positionnement approprié à l'œil.

Le Microscope Zoom Stéréo Binoculaire d'Olympus SZ61 dispose d'un rapport de zoom de 6,7:1 avec un grossissement variable de 6,7 - 45X lorsqu'il est utilisé avec un oculaire de 10X Le Microscope Zoom Stéréo Trinoculaire SZ-61R est également disponible et dispose d'un tube trinoculaire avec un objectif de monture C 0,5X intégré, facilement compatible avec des [caméras](#) numériques ou vidéo.

Informations techniques

Description		Stock No.
Microscope Body (One Required)		
SZ51 Stereo Microscope Body with 45° Binocular		#88-124
SZ61 Stereo Microscope Body with 45° Binocular		#88-125
SZ61 Stereo Microscope Body with 45° Trinocular (C-Mount)		#88-126
Stands (One Required)		
Basic Reflected Light Stand with Reversible Stage Plate		#88-207
Transmitted and Reflected LED Light Stand		#88-208
Eyepieces (Two Required)		
10X Widefield Eyepieces		#88-108
15X Widefield Eyepieces		#88-109
20X Widefield Eyepieces		#88-110
30X Widefield Eyepieces		#88-111
Objectives (One required)		
0.3X Objective Lens		#88-127
0.4X Objective Lens		#88-128
0.75X Objective Lens		#88-129
1.5X Objective Lens		#88-130
2X Objective Lens		#88-131
Optional Accessories		
Ball & Socket Stage for Basic Light Stand		#88-123

Microscope Specifications									
Microscope Body	Zoom Magnification	10XEyepiece #88-108		15XEyepiece #88-109		20XEyepiece #88-110		30XEyepiece #88-111	
		F.N. 22		F.N. 16		F.N. 12.5		F.N. 7	
		Mag.	FOV	Mag.	FOV	Mag.	FOV	Mag.	FOV
SZ61	0.67X	6.7X	32.8mm	10.1X	23.9mm	13.4X	18.7mm	20.1X	10.4mm
	1X	10X	22mm	15X	16mm	20X	12.5mm	30X	7mm
	2X	20X	11mm	30X	8mm	40X	6.3mm	60X	3.5mm
	3X	30X	7.3mm	45X	5.3mm	60X	4.2mm	90X	2.3mm
	4.5X	45X	4.9mm	67.5X	3.6mm	90X	2.8mm	135X	1.6mm
SZ51	0.85X	8X	27.5mm	12X	20mm	16X	15.6mm	24X	8.8mm
	1X	10X	22mm	15X	16mm	20X	12.5mm	30X	7mm
	2X	20X	11mm	30X	8mm	40X	6.3mm	60X	3.5mm
	3X	30X	7.3mm	45X	5.3mm	60X	4.2mm	90X	2.3mm
	4X	40X	5.5mm	60X	4mm	80X	3.1mm	120X	1.8mm

Olympus SZ51/SZ62	Olympus SZX7
Affordably priced with high resolving power	Larger zoom ratio than the SZ51/SZ61
Utilizes a Greenough optical system	Brighter, higher resolution images
	Utilizes an advanced Galilean optical system
